

## Next Generation Sequencing (NGS)

Clonal sequencing using the Illumina MiSeq platform has been introduced as a diagnostic service in the Laboratory in June 2013 replacing conventional Sanger sequencing of Glycogen Storage Diseases as required following discussion with referring clinicians. Our validation studies\* demonstrated a test sensitivity and specificity comparable to conventional sequencing. The laboratory will be introducing a number of additional services utilising this technology. In addition, the laboratory will be introducing a number of clonal sequencing services using the LifeTechnologies Ion PGM platform.

\* In-house data

## Current NGS Services

### Glycogen Storage Disorders

(Design ID: 0423281).

<b>Primary liver panel</b>		<b>Liver panel</b>	
<b>(available during validation phase)</b>	AGL		AGL
	FBP1		FBP1
	G6PC		G6PC
	GBE1		GBE1
	GYS2		GYS2
	PHKA2		PHKA2
	PHKB		PHKB
	PHKG2		PHKG2
	PYGL		PYGL
	SLC37A4		SLC37A4
			GYG2
			PFKL
			SLC2A2
<b>Primary muscle panel</b>		<b>Muscle panel</b>	
<b>(available during validation phase)</b>	AGL		AGL
	GAA		GAA
	GBE1		GBE1
	GYS1		GYS1
	PFKM		PFKM
	PGAM2		PGAM2
	PHKA1		PHKA1
	PHKG1		PHKG1
	PYGM		PYGM
			ALDOA
			ENO3
			FBP2
			GYG1
			LDHA
			PGK1

# Testing Workflow

## Library Preparation

- Shearing of genomic DNA using the Covaris E220 sonicator.
- End repair, A tailing and ligation of adaptors using SureSelectXT library system (Agilent Technologies).
- Enrichment by SureSelect target enrichment (Agilent Technologies) using custom in house designed probes. Samples have barcode tags added following target enrichment.
- Quantification of tagged samples by QPCR using KAPA Library quantification kit for Illumina and lightcycler 480.
- Sequencing on the Illumina MiSeq using the MiSeq Reagent Kit v2 performing 2 x 150bp end paired reads.

## Data Analysis

- Data analysis using NextGENe v2.3.1 & v2.3.3 (SoftGenetics) software.
- Data aligned to human genome build hg19 (dbSNP135).
- Generation of mutation reports and coverage curves.
- Read depth checked visually using Alamut v 2.2 (Rev1) (Interactive Biosoftware)
- Variants filtered following Best Practice Guidelines for the evaluation of pathogenicity and the reporting of sequence variants in clinical molecular genetics (Association for Clinical Genetic Science).
- Variant frequency checked using dbSNP build 137.
- A minimum threshold of 30-fold read depth is set for diagnostic sequencing.

## Post analysis

- Confirmation of all clinically significant sequence variants by Sanger sequencing.
- Creation of a diagnostic report combining clonal and Sanger sequence data.